

Table S1. Elution method via HPLC-DAD-FLD for the quantification of CBD and NaF.

Time (min)	Acetonitrile (%)	Acid Water (pH 3.2)
0.1	20	80
2	25	75
10	50	50
11	50	50
20	99	1
35	20	80

Table S2. Elution method via HPLC-MS for the quantification of CBD.

Time (min)	Acetonitrile (%)	Acid Water (pH 3.2)
0.1	25	75
0.5	25	75
8	100	0
14	100	0
15	25	75
20	25	75